

submit a quote by visiting our website
or calling toll free 888.551.4062
mention **PROMO CODE: summer10**
to receive a discount
www.shor-line.com



Single condo configurations starting at \$486.



Laminate shown above.

CONDOS

{long service life & healthier environment}

**BEAT THE
\$ PRICE \$
INCREASE**

Stainless Steel shown below.

Stainless Steel Cat Suite
Multiple Patents Pending



902.4830.20 Single Unit
902.0103.23 Double Unit with Mobile Platform

Benefits:
Stainless steel construction provides the best surface for disinfection and an incredibly long service life. All features are designed to relieve the stress of the occupant, thus leading to a healthier environment and safe habitat for both feline and handler.

Features:
• Stainless steel door design discourages paw entrapment.
• Lockable door latch to prevent escapes.
• Silent Latch Mechanism eliminates noise and creates a stress-free environment for patients.
• Integrated compartment divider can be removed for cleaning.
• Portals can be opened and closed without opening cage doors or disturbing occupants.
• Litter compartment barrier keeps litter from falling onto the floor.



{902.0103.23}



{Silent Latch Mechanism}



{Open & Close Portals}

GIVEAWAY

{your pick: one Quiet-Time Kennel Door Cover or Card Holder*}

*With any single cage purchase receive your choice of one Quiet-Time Kennel Door Cover or Card Holder (card holder shown on accessory page). While supplies last.



{179.2424.00}



Quiet-Time Kennel Door Covers

179.2424.00 24" w x 24" h kennel door
179.2430.00 24" w x 30" h kennel door
179.3030.00 30" w x 30" h kennel door
179.3630.00 36" w x 30" h kennel door

Benefits:

Are you hanging old towels over your kennel doors to reduce stress and separation anxiety? Most cats and dogs become calm when you block their visibility of people and animals. Quiet-Time Kennel Door Covers also provide some isolation and solitude for patients during treatment or recovery.

Features:

- Color-fast blue nylon construction is stain and snag resistant.
- 3" deep pocket fits snugly over the top of door and secured with plastic snaps.
- Clear vinyl window to display client/patient information.
- Machine washable.

Available for purchase.

Like us on Facebook.
 Follow us on Twitter.

Promotional discount is not limited to products shown in this brochure to products listed in this brochure. All products listed in this brochure will be invoiced at the discount price communicated during the promotion period. Prices are F.O.B. Kansas City, KS. Brochure promotional discount cannot be used in conjunction with any other promotional discount on Shor-Line products. Prices are subject to change without notice. Prices do not include shipping & handling charges or any applicable taxes. Shor-Line is not responsible for any typographical errors. Complete Shor-Line terms and conditions are available on our website: www.shor-line.com. Promotional discount valid in US only. Sale effective now thru 7/30/10.



PROMO CODE: summer10

call 888.551.4062

fax 913.281.5339

to place an order

www.shor-line.com

Shor-Line/Schroer Manufacturing Co.
511 Osage Avenue
Kansas City, KS 66105 U.S.A.

SHOR-LINE

SUMMER PROMOTION

{now thru July 30, 2010}



{Fold-Up Wall Mount Exam Table Attached to Stainless Steel Chase}

**BEAT THE
\$ PRICE \$
INCREASE**

Promotional discount is not limited to products shown in this brochure. Visit www.shor-line.com for a complete product listing.

Cages, Condos, Runs, Grooming, Cabinetry & Accessories

Outfit your facility with top of the line products offered only by Shor-Line. We offer a complimentary design service, financing options and a committed representative to answer all your animal care equipment questions. Mention **PROMO CODE: summer10** and receive a discount on your next animal care equipment order!



SHOR-LINE

{Unleashed Innovation}



Join our community. With Shor-Line's My Account feature you will be able to create a wish list, obtain quotes, review past quotes and have access to valuable resources. Set up your profile today! Visit www.shor-line.com.

PRSRT STD
U.S. POSTAGE
PAID
KANSAS CITY, MO
PERMIT NO. 6744

**BEAT THE
\$ PRICE \$
INCREASE**

CAGES

TYPE
SHOR-LINE CAGES

{dare to compare}

We have been manufacturing cages for decades. We have customers who's cages have outlasted multiple practice owners. Our products will out shine any competitor based on quality, craftsmanship and materials. We dare you to compare.

Assemblies starting as low as \$1388.



Cage Feature	Shor-Line Type 3 Cages	Competitor A	Competitor B
Door			
Latch			
Hardware			
Hinge			
Frame			

Trim Hardware

Solid stainless steel with rounded corners-no sharp edges or points.

Patented Latch

Designed to be durable, secure and quiet. Acetal bushings provide smooth operation and sound dampening qualities.

Sound Dampening

New sound and vibration absorbing technology. Standard equipment on all Shor-Line cages.

World's Mightiest Cage Door-Bar None...

The one piece frame and wrap around hinge/latch bars are formed from 3/8" diameter solid stainless steel rods. All additional horizontal rods are 1/4" in diameter and all vertical rods are 3/16" in diameter. Every intersection of horizontal and vertical rods are welded, providing incredible strength. Electro polishing is the final process resulting in a smooth and bright finish that looks great and is easy to clean...and practical too! Loops on the door prevent paper, pads or towels from sliding out the front of the cage.

Fully Rounded Corners & Fluid Retaining Lip

A smooth rounded cove transitions the back and all internal sides of the cage. A slightly raised, fluid retaining lip prevents drainage out the front of the cage.

New Stainless Steel Tube Frame

Front door frames are heavy duty 1" square stainless steel tube to provide you many years of solid structural support for your cage assemblies. All edges and corners on our frames have been radiused and de-burred.

Patent Pending Hinge Design

The hinge is constructed of 14 gauge austenitic stainless steel that is encapsulated by an injection molded, fiberglass reinforced nylon 6/6 thermoplastic polymer.

submit a quote by visiting our website or calling toll free 888.551.4062 mention **PROMO CODE: summer10** to receive a discount www.shor-line.com

RUNS

{complimentary design service}



Run Gates starting at \$399. Build your assembly from items below!

{great deal on overstock kennels*}

*Listed below.

971.4872.01	4'w x 6'h	Stainless Steel Gate
961.7172.42	6'1 x 6'h	Solid Side Panel
961.4772.40	4'1 x 6'h	Solid Side Panel
961.7172.49	6'1 x 6'h	One Piece Galvanized Side Panel
961.4673.56	4'w x 6', 1.5'h	Back Panel One Side PVC
961.4642.91	4'w x 4'1	Double Decker Upper Floor
941.4645.36	4'w x 4'1	Double Decker Wire Mesh Top
961.4673.53	4'w x 6', 1" h	Double Decker Back Panel One Side



888.551.4062

GROOMING

{top notch styling with affordable bathing products}

Step-Up Bathing Tub with Ramp/Steps



{Steps Feature}

Easily transitions from stairs to a ramp in no time!



{Ramp Feature}



{Anti-Fatigue Mat}

Bathing and grooming animals demands a dedicated space and professional equipment. Visit www.shor-line.com to view more Shor-Line grooming products.

CABINETS

{quality craftsmanship & custom design}



BEAT THE \$ PRICE \$ INCREASE

Cabinetry adds a residential comfort that is a welcoming change from the clinical atmosphere. Patients' owners can rest easy knowing that their pet is being cared for at an extension of their home. Call today to speak with a committed representative about your next project!

888.551.4062

ACCESSORIES

{complete your animal care area with these must haves}

Card Holders



802.0000.01 3" x 5" Stainless Steel \$8

Patient Stretcher



918.0003.04 \$130

Animal Handling Gloves



918.0001.01 8" Gauntlet Large \$55
918.0001.02 8" Gauntlet Small \$55
918.0001.03 19" Gauntlet Large \$94
918.0001.04 19" Gauntlet Small \$94

Woods/Exam Hand-Held Light with Magnifier



913.2004.00 \$349

Feline Scale



905.0400.01 \$279

Versatile Bird Perch



905.0400.03 \$20

Economy Lab Stool



918.0004.34 \$94

Small Animal/Avian Scale



905.0400.02 \$169