

BLUE-LINE FOLD UP WALL MOUNT TABLES

Fold-Up Wall Mount Exam Table

Benefits:

Maximize floor space and style with the Blue-Line Fold-Up Wall Mount. It provides both a sturdy exam space and ample surface area, all without getting in the way of your day-to-day exam activities.

Features:

- The ONLY Fold-Up Wall Mount with radiused corners.
- Laminated polyurethane edging holds up to one liter of fluid.
- Torsion spring technology allows you to easily and safely raise and lower the table
- Weight Capacity: 200 lbs.
- Cabinet mounting options available for extra storage.
- Specs: The table top is 44½”L x 22½”W x 1”H, and it sits five inches from the wall when folded up. Easily mounts to any standard 2” x 4” reinforced wooden stud wall with 16” or 24” stud centers.
- Compatible with our Stainless Steel Utility Chase.



- 903.1530.04 Lateral, Standard Blue
- 903.1530.05 Lateral, Jade
- 903.1530.06 Lateral, Taupe Gray
- 903.1530.07 Lateral, Black
- 903.1531.04 Peninsula, Standard Blue
- 903.1531.05 Peninsula, Jade
- 903.1531.06 Peninsula, Taupe Gray
- 903.1531.07 Peninsula, Black

Additional Colors:

Standard Blue Black Jade Taupe Gray

