

FELINE HOUSING

R E I M A G I N E D



feline comfort SUITE

INTRODUCING A NEW VISION FOR FELINE HOUSING

We are pleased to announce our new cat condos that are strong enough to withstand the harshest cleaning liquids and still have a silky, smooth surface for cat paws.



STRENGTH THROUGH TECHNOLOGY

We craft the condos from a revolutionary fine-celled foam board created with a new manufacturing process. This produces a hard but smooth outer structure that has a strong cellular core.

What that means to your feline clients is stronger than wood and impervious to water and urine. Its strength can take a pounding from the largest cats without damage. The silky matte finish has just the slightest texture for secure footing.



We reached out to the construction industry to use the latest technology to create durable hinge that will stand up to repeated opening and closing. It is an injection-molded, fiberglass-reinforced nylon thermoplastic polymer.

LUXURY FELINE VETERINARY BOARDING

We've tried to get inside the feline brain through talks with our customers and thought leaders as well as time spent observing cats. We purposely reimagined how to make cats comfortable when spending time away from home in a kennel. The result is the Feline Comfort Suite that is easy on your eyes and soft on their paws.

Two door options are available:
The **clear, polycarbonate door** covers the entrance while providing diagonal ventilation—a shape that makes cleaning easier. Three adjustable polymer vents in the back keep air flowing through the unit. All provides a clear view for cat and care takers.

The **Tri-Panel** combines a grill front with a polycarbonate band in the middle that gives the cat a bar-free view at eye level. The grill gives good ventilation through the unit.

The unit is crafted from durable, fine-celled foam board that is stronger than wood wood and impervious to moisture.

Large rods allow ventilation in the litter area. The unique divider design has a cleaner edge that removes possible debris from the slot. You can adjust "open" or "closed" portal status from outside the suite.

Litter and fluid retention strip at the cage front helps avoid drips.

The state-of-the-art injection-molded, fiberglass-reinforced nylon thermoplastic polymer hinge offers durability in a quiet, lightweight product.

We crafted the latch to inset into the cage side panel to protect gentle paws and increase latch life. The Stainless Steel latch bar is lockable.

Litter space gives ample vertical room for a cat to take care of business in comfort. The polycarbonate section on the bottom has a clip for attaching records or information. The grill provides good air circulation.



U.S. Patent #7,805,808



LITTER AREA (One Size)

Inside dimension: 16½"W x 28½"H x 28"D
Outside dimension: 18"W x 30"H x 28¾"D

LIVING AREA

(Two Sizes Available)

Medium

Large

Inside dimension:	22½"W x 28½"H x 28"D	28½"W x 28½"H x 28"D
Outside dimension:	24"W x 30"H x 28¾"D	30"W x 30"H x 28¾"D

A NEW ERA OF FELINE BOARDING COMFORT HAS BEGUN

Our Feline Comfort Suite was developed with the latest feline housing and low stress design principles. Notice the many premium features that make the product more appealing to cats as well as easier for caretakers.

We are pairing our craftsmanship with a revolutionary fine-celled foam board created with a new manufacturing process. The result is a feline housing unit that will become as legendary and durable as our Stainless Steel cages.



ACCESSORIZE

Optional Kat Kave: Two-level cave lets cats lounge on top, snuggle in the bottom or wind through the three openings below. Made of durable, easy-to-clean black, high molecular weight polyethylene. SPECS: 21½"W x 12"H x 12¼"D.

Privacy Quiet-Time Cover. Colorfast nylon cover snaps onto the cage to provide a space of privacy for the cat. It comes in blue and is machine washable. Note: Designed for Feline Comfort Suites with grill fronts. Does not fit polycarbonate door. SPECS: 12"W x 20½"H



800.444.1579

SHOR-LINE.COM

SHOR-LINE®